



## How much blood is too much to draw from a neonate?

In these tiny patients, any amount seems like too much. But when blood is needed, **the least amount possible** is the goal.

### CLINITUBES and safeCLINITUBES

Accurate, clot-free testing for pH, blood gas, electrolytes, oximetry and metabolites on 65  $\mu$ L

Radiometer's CLINITUBES and safeCLINITUBES™ capillary samplers are ideal for capturing small-volume samples from neonates and infants. safeCLINITUBES are unbreakable, making them OSHA compliant and safe for both baby and clinician in the NICU.

Available in balanced or high heparin concentrations, preheparinized CLINITUBES and safeCLINITUBES prevent biased electrolyte values and ensure accurate, clot-free testing for pH, blood gases, oximetry and metabolites.



150  $\mu$ L



100  $\mu$ L



Just 65  $\mu$ L with CLINITUBES



[www.radiometeramerica.com/clinitubes](http://www.radiometeramerica.com/clinitubes)

Unbreakable plastic • OSHA-compliant • Pre-heparinized • Clot catchers available • Ideal for neonates and infants